

=====

Ref. GCDMS103

=====

## Description

=====

## Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE            | Medical device | Class I    |
| Bar code          | 3700468208136 | HS Code        | 9018908400 |
| Warranty (months) | 24            | Country        | France     |

### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES