

PATIENT STOOL WITH NYLON BASE 5 DOUBLE CASTORS BLUE

Ref. OS007_L



PATIENT STOOL WITH NYLON BASE 5 DOUBLE CASTORS BLUE



Ref. OS007_L

Plastic base stool (seat 3.5 cm)

Description

- Seat ø: 43 cm, padded in very high density foam (30 kg / m3)
-
- Fireproof seat (fire resistance classification M1)
- Axis, tube ø: 5.5 cm in stainless steel
- Black plastic base with 5 branches ø: 58 cm, very stable with double casters
- Height adjustable by simple pressure and gas piston
- Height variation from 53 to 70 cm
- Maximum load: 150kg

Matching stool to the OSA003 sofa (several colors available for both models)

- Available in Black (OS007_N) and Lavender (ref: OS007_L)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040196138	HS Code	9402900000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	17	Packed unit / Height (cm)	42
Packed unit / Length (cm)	41	Packed unit / Weight	6.2
Bare Unit / Weight (kg)	5.2		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Lavender	Packed unit / Quantity	1
--------	----------	------------------------	---