

Ref. 522711



30/09/2026 01:23:56 – For information purposes only, can be subject to change

=====

Ref. 522711



Description

Cuff available in :

Adult (M), nylon, velcro, black (522711)

Adult (L), nylon, black (522712)

Child/small adult (S), nylon, black (522713)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028176	HS Code	9018199560
Warranty (months)	60	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.092

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---