


Ref. 337843



06/10/2026 00:00:58 – For information purposes only, can be subject to change

=====

Ref. 337843



Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Spengler hand-held BPMs Protocol patient monitors

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446026943	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.125
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 06/10/2026 00:00:58 – For information purposes only, can be subject to change