


=====

Ref. 337845



 05/10/2026 22:33:10 – For information purposes only, can be subject to change

=====

Ref. 337845



Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Spengler hand-held BPMs Protocol patient monitors

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028404	HS Code	9018901000
Warranty (months)	6	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.1

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---