

LUXASCOPE OPHTHALMOSCOPE LED 3.7 V - LUXAMED®

Ref. F1.426.114

LUXAMED®



Rechargeable USB



07/10/2026 01:34:10 – For information purposes only, can be subject to change

LUXASCOPE OPHTHALMOSCOPE LED 3.7 V - LUXAMED®

LUXAMED®Ref. F1.426.114

Direct ophthalmoscope in pocket size

Description

Direct ophthalmoscope in pocket size

Correction values: + 20 bis - 20 D ($\pm$ 1 2 3 4 6 8 10 15 20 D)

5 different apertures: reticle, large circle, small circle, semicircle, green filter

Illumination: LED illumination

Can be dimmed for an individual glare-free brightness

Automatic switchoff after 3 minutes

Life expectancy of LEDs: approx. 50.000 hours by constant current drive

Colour temperature: approx. 4.200 K

Optimized design of head

Two-component handle (aluminium/fibre-glass reinforced), also combinable with LuxaScope LED otoscope- and dermatoscope head

Electronic ON/OFF push-button switch

Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive

Incl. recharger with USB plug (adapter EU/UK/US)

Available in :

- Black: F1.426.114

- White : F1.426.914

LUXASCOPE OPHTHALMOSCOPE LED 3.7 V - LUXAMED®

LUXAMED®Ref. F1.426.114

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964339	HS Code	90185090
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Black
--------	-------